

Diagnostic Report Format & Templates

Diagnostics, Labs & Imaging

Document type: Plan / SOP / strategy / guide · Generated 14 Oct 2026

This guide explains what the tool produces, how to use the output, and what to watch for. The actual output (draft/report) is generated on the tool page.

Background

Designs diagnostic report formats and templates — layout, patient and test details, reference ranges, flags, interpretation, signatory and disclaimers. Use it to standardise reports. For critical results use the Critical Value Alert Protocol.

Purpose of this tool & output

Design clear, compliant report formats

How the output is used

Use it to design clear, compliant report formats — templates with reference ranges, flags and the signatory details a report needs.

Important cautions

Reports must carry correct patient identifiers, method-appropriate reference ranges, and a genuinely qualified signatory block. Include limitations and advise clinical correlation. Align formats with NABL and PC-PNDT (for imaging) requirements.

Companion documents

- Test Menu & Catalogue Design
- Health Checkup Package Designer
- Pre-Analytical Sample Workflow SOP
- Phlebotomy & Sample Collection SOP
- Turnaround Time (TAT) Management
- Outsourced / Referral Test Plan

This is AI-assisted guidance. Verify with a qualified professional before any binding decision. Fees, authorities and processes vary by state and change over time — confirm with your state authority.